



SURU INTERNATIONAL PVT. LTD.

MONOFILAMENT POLYPROPYLENE SUTURE

Product Name: Monofilament Polypropylene Suture - SURULENE®


Item Code: 1013

Features :

- Product : Monofilament Polypropylene Suture SURULENE® (Non-Absorbable Suture U.S.P)
- Material : Monofilament Polypropylene
- Colour : Blue
- Dye Used : CU-Phtalocyanine blue (US FDA approved dye)
- Coating : The Product does not require any coating and is sold untreated

Available Configuration :

Size (U.S.P)	6-0	5-0	4-0	3-0	2-0	1-0 or 0	1	2
Metric	0.7	1	1.5	2	3	3.5	4	5



NON - ABSORBABLE SURGICAL SUTURE U.S.P. **P1A70DZ30**
Monofilament Polypropylene Suture

1 (4 Metric)
70 cm














1/2 Circle
 Round Bodied
 30mm

STERILE EO 35°C





















MFG. LICENCE # KD 459
LOT : L18015
AWT : 2018.10
EXP : 2023.09


















Manufactured by :
SURU INTERNATIONAL PVT. LTD.
 Plot # 5, Poonamia Comp., Saravali, Jawhar Rd.,
 Dahanu 401 602, MAH., INDIA • www.suru.com
 Customer Care No.: + 91 22 2820 0707. Email : suru@suru.com

SURULENE®

NEEDLE LENGTH	NEEDLE DESCRIPTION	SUTURE LENGTH	SUTURE SIZE (USP)				
			10-0 (0.2)	9-0 (0.3)	8-0 (0.4)	7-0 (0.5)	6-0 (0.7)
 4 mm	3/8 Circle MP Spatula	45 cm	SU A869				
 4 mm Double Armed	3/8 Circle MP Spatula	45 cm	SU A869D				
 6 mm	3/8 Circle MP Spatula	45 cm	SU A868	SU A875	SU A871		
 6 mm Double Armed	3/8 Circle MP Spatula	45 cm	SU A868D	SU A875D	SU A871D		
 8 mm	3/8 Circle MP Round Bodied	45 cm				SU A873	
 8 mm Double Armed	3/8 Circle MP Round Bodied	45 cm				SU A873D	
 8 mm Double Armed	3/8 Circle Round Bodied	70 cm					SU 893
 9 mm Double Armed	3/8 Circle Taper Cutting	70 cm					SU 874
 10 mm Double Armed	3/8 Circle Round Bodied	70 cm					SU 895
 10 mm	3/8 Circle Cutting	70 cm					SU A817
 10 mm Double Armed	3/8 Circle Cutting	70 cm					SU A817D
 16 mm	Straight MP Spatula	45 cm	SU AY868				
 16 mm Double Armed	Straight MP Spatula	45 cm	SU AY868D				

Abbreviation : MP : Micro Point

NEEDLE LENGTH	NEEDLE DESCRIPTION	SUTURE LENGTH	SUTURE SIZE (USP)				
			6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
 12 mm Double Armed	3/8 Circle Round Bodied	70 cm		SU 827			
 12 mm	3/8 Circle Reverse Cutting	45 cm	SU 878		SU 879		
		70 cm		SU 880			
 13 mm Double Armed	3/8 Circle Round Bodied	70 cm	SU 829				
 15 mm	3/8 Circle Slim Blade Cutting	70 cm	SU 823	SU 826			
 16 mm	3/8 Circle Round Bodied	70 cm		SU 881			
 16 mm Double Armed	3/8 Circle Round Bodied	70 cm		SU 882			
 16 mm	3/8 Circle Cutting	70 cm			SU 870		
 16 mm Double Armed	1/2 Circle Round Bodied	90 cm			SU 850		
 16 mm	1/2 Circle Round Bodied	70 cm			SU 849		
 17 mm	1/2 Circle Taper Cutting	90 cm			SU 857		
		70 cm			SU 899		
 17 mm Double Armed	1/2 Circle Taper Cutting	90 cm			SU 853	SU 836	SU 837
		70 cm			SU 835	SU 897	SU 898
 19 mm	3/8 Circle Cutting	75 cm			SU 834		
 20 mm	1/2 Circle Round Bodied	70 cm				SU 887	
 20 mm	3/8 Circle Cutting	45 cm			SU 894	SU 893	
 20 mm Double Armed	1/2 Circle Round Bodied	90 cm			SU 861		SU 861A
 20 mm Double Armed	1/2 Circle Cutting	90 cm					SU 861B
 22 mm	3/8 Circle Cutting	70 cm				SU 018	SU 017
 25 mm	1/2 Circle Round Bodied	90 cm				SU 825	
		70 cm				SU 838	
 25 mm Double Armed	1/2 Circle Round Bodied	90 cm				SU 852	SU 852A
		70 cm				SU 822	
 25 mm Double Armed	1/2 Circle Cutting	90 cm					SU 855

NEEDLE LENGTH	NEEDLE DESCRIPTION	SUTURE LENGTH	SUTURE SIZE (USP)				
			4-0 (1.5)	3-0 (2)	2-0 (3)	0 (3.5)	1 (4)
	3/8 Circle Cutting	45 cm	SU 856	SU 856A			
	1/2 Circle Taper Cutting	90 cm			SU 806		
		70 cm			SU 807	SU 805	
	1/2 Circle Taper Cutting	90 cm	SU 867		SU 877		
		70 cm			SU 808		
	1/2 Circle Round Bodied	90 cm				SU 824	
	3/8 Circle Cutting	45 cm		SU 858			
	1/2 Circle Round Bodied	90 cm			SU 885		
	1/2 Circle Round Bodied	70 cm			SU 841	SU 842	
	1/2 Circle Round Bodied	70 cm				SU 830	SU 840
	1/2 Circle Taper Cutting	90 cm				SU 891	
		100 cm					SU 896
	1/2 Circle Round Bodied	100 cm				SU 886	
	1/2 Circle Cutting	100 cm				SU 895	
	1/2 Circle Round Bodied	70 cm					SU 843
	1/2 Circle Round Bodied	90 cm				SU 845	
		100 cm					SU 892
	3/8 Circle Cutting	100 cm					SU 890
	1/2 Circle Reverse Cutting	100 cm					SU 883
	1/2 Circle Round Bodied	150 cm (loop)					SU 834
	Straight Cutting	70 cm		SU 888			

<http://www.suru.com/Surgical-Sutures/SURULENE-Monofilament-Polypropylene-Suture>